vyvo™


 **sense**

The Smartband that
heightens your Senses,
listening to your body's voice.

QUICK START GUIDE

EN

1. PACKING LIST


The box contains the following items:

SENSE

Wristband

Instructions



2. HOW TO OPERATE

2.1. Charging

As shown below, press the release button on the back of the bracelet to make the main device detach from the strap.

Plug the detached main device into the USB port for charging (only one end with the metal electrode can be charged). When the device starts charging, there will be a vibration prompt, and the screen will show the status of charging.

Input voltage: 5V DC.

Charging time: about 1.5 hours.

Note: Please use a certified regular charger when charging.

2.2. Power On

Press and hold the screen for power on. If the device does not power on, please fully charge the device first.

2.3. APP Download

Download the VYVO Smart App from the App Store or Google Play Store.

Note: To use "VYVO Smart App" on an iPad, you need to select "iPhone only" when downloading it from the App Store.

2.4. Compatibility & Support

Support Android 4.4 or above, iOS 9.0 or above.

2.5. Connection

After the download has finished, open the App and register an account, then complete the connection process according to the operating instructions on the App.

NOTE:

a. In order to successfully connect a phone to the smartband, you

need to turn on the Bluetooth on the phone and connect to your device through the App.


b. The first time you connect to the App, the time on the smartband will not be consistent with your actual local time. Once the connection has been completed the time along with all the data for steps, calories & distance on the smartband will be reset.

c. The app needs permission to access your location if you are using an Android phone, so you need to ALLOW LOCATION ACCESS in System Settings in order to connect successfully.

2.6. Wearing the Device

To make the measured values more accurate, please wear the device one finger away from your wrist bone.

When exercising, please wear close to your skin and make sure that the device will not move along, slide, or jiggle with your movement.

3. INTRODUCTION FOR IMPORTANT FEATURES

3.1. Heart Rate Monitoring

The device contains PPG heart rate sensors and combined with HR algorithm, accurately monitors your heart.

Five heart rate zones are displayed during the exercises, but the detailed data can be checked after connecting and syncing to the VYVO Smart App.

Note: The light signal transmission might be blocked if your skin is too dry or too dark or with excess hair, which may lead to a measurement failure.

3.2. 15 Sports Modes

Slide up or down from the home screen to access "Sports" and there are 15 types of exercise for your choice: running, cycling, walking, treadmill, basketball, badminton, table tennis, tennis, football, volleyball, jumping rope, sit-ups, plank, aerobics, and swimming. After the exercise, you can directly check your results on the smartband. Detailed data shall be viewed after connecting and syncing to the app.

3.3. Sleep Monitoring

Wear the smartband when you go to bed at night, then you can check your sleep data monitored on the app after you get up in the morning. The smartband goes into sleep monitoring from 8:00 pm to 9:00 am the next day.

Note:

- 1.The sleep monitoring function will be stopped after you get up and move for 5-10 minutes.
2. Daytime sleep data won't be recorded.

3.4. Change the Smartband Face

Slide up or down from the home screen to access "Set", click "Set" and find "Dial" (Press and hold it to find more dials, swipe up or down to switch dial and click the bottom of the screen to confirm).

3.5. Alarm/Schedule

Alarm clock or Schedule are all set within the app, and the setting will be automatically synchronized to the device.

4. MORE INFORMATION

4.1. Waterproofing instructions

The waterproof level of SENSE is IP68.

Note:

1. Swim, shower, and wash hands only under cold water.
2. The device does not have a waterproof effect on corrosive liquids such as seawater, acid & alkaline solutions, and chemical reagents. In case of accidental contact with corrosive liquid, please clean it with clear water and dry it. Damages or defects caused by abuse or improper use are not covered by the warranty.

4.2. Firmware Upgrade

Go to the Device interface of the app and select firmware upgrade.

Notes:

1. Make sure that the battery level is more than 50% before upgrading.
2. During the upgrade process, you cannot quit halfway if the progress bar moves, keep the screen of your phone bright, and only when the upgrade is completed can you exit the interface. Otherwise, the upgrade will fail. Wait for the device to restart automatically if the upgrade fails. Then reconnect to the app for upgrading again.

4.3. Specification parameter

Product model	SENSE
Size	37 × 20 × 12 mm
Battery Capacity	Li-Polymer battery 75 mAh
Battery life	5-7 days (Exercise an average of half an hour a day; Receive an average of 30 messages and 4 calls a day; Raise your wrist to make the screen turn on 50 times; 24h heart rate monitoring turned on).
Working condition	0-40°C
Weight	About 27g
Display Size	0.96in TFT
Screen resolution	80 x 160 Pixels

5. DEVICE MAINTENANCE

5.1. Device care

Notice:

1. Do not use a sharp object to clean the device.
2. Avoid using solvents, chemical cleaners, or insect repellents that might damage plastic components of the device.
3. Thoroughly rinse the device with fresh water after exposure to chlorine, salt water, sunscreen, cosmetics, alcohol, or other harsh chemicals. Prolonged exposure to these substances can damage the device.
4. Avoid extreme shock and harsh treatment since it can reduce product life.
5. Do not expose to extreme temperatures that are too cold or too hot for a long time, which can cause permanent damage.
6. After each training, please rinse the smartband with clean water.

5.2. Cleaning the Device

1. Gently wipe the device by using a cotton cloth with a mild neutral detergent.
2. Wait for it to dry.


Note: Even the fainter sweat or moisture may cause the corrosion of the charging terminal when charging the device, which will also hinder the transmission of data and affect the charging.

6. IMPORTANT SAFETY AND PRODUCT INFORMATION

- If you have a pacemaker or other internal electronic device, consult your physician before using a heart rate monitor.
- The inside optical wrist heart rate monitor will occasionally emit green light and flashes. Consult your physician if you have epilepsy or are sensitive to flashes of light.
- Always consult your physician before beginning or modifying any exercise program.
- The device, accessories, heart rate monitor, and related data are intended to be used only for exercise monitoring, not medical purposes.
- The heart rate readings are for reference only, and no liability is accepted for the consequences of any wrong interpretation.
- Although the optical wrist heart rate monitor technology typically provides an accurate data for user's heart rate, there are inherent limitations with the technology that may cause some of the heart rate readings to be inaccurate under certain circumstances, including the user's physical characteristics, installation of the device and type & intensity of activity.

- The activity trackers rely on sensors to track your movement and other metrics. The data and information provided by the device are intended to assess your activity and metrics closely, but may not be completely accurate, including the data for step, sleep, distance, heart rate, and calories.
- Do not immerse or expose batteries to water or other liquids.
- Do not expose the device to a heat source or in a high-temperature location, for example, in an unattended vehicle in the sun. To prevent the possibility of damage, remove the device from the vehicle or store it out of direct sunlight, such as the glove box.
- When storing the device for an extended period of time, please store within the temperature ranges specified in the printed manual, which included in the product packaging.
- Do not use a power and/or USB cable that is not approved or supplied by the manufacturer.

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